

IHE Integration Statement

Vendor	Product Name	Version	Date
FUJIFILM Corporation	FDX Console / Console Advance (DR-ID 300CL)	8.1	25/12/2014
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actor Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist Query Modality Group Case	
Patient Information Reconciliation (PIR)	Acquisition Modality		
Radiation Exposure Monitoring (REM)	Acquisition Modality		
Portable Data for Imaging (PD)	Portable Media Creator		
Internet address for vendor's IHE information : http://www.fujimed.com/sub/product_bulletins.html			
Links to Standards Conformance Statements for the Implementation			
HL7	-		
DICOM	http://www.fujimed.com/sub/product_bulletins.html		
Links to general information on IHE			
http://www.jira-net.or.jp/ihe-j (Japan) http://www.ihe.net (North America) http://www.ihe-europe.net (Europe)			