

SYNAPSE

PACS

Synapse PACS IHE Integration Statement

- * SYNAPSE is a registered trademark of FUJIFILM MEDICAL SYSTEMS, U.S.A., Inc. in United States, Australia and Canada.
- * SYNAPSE is a registered trademark of FUJIFILM Corporation in China, Japan, Austria, Belgium, Denmark, Finland, France, Germany, Greece, Ireland, Italy, Luxembourg, The Netherlands, Portugal, Spain, Sweden, United Kingdom, Latvia, Slovakia, Slovenia, Cyprus, Lithuania, Hungary, Czech Republic, Malta, Poland, Estonia, Bulgaria, & Romania.

Manufactured by:

FUJIFILM Medical Systems U.S.A., Inc.
3020 Carrington Mill Blvd, Suite 500
Morrisville, NC 27560 U.S.A.

Manufactured for:

FUJIFILM Medical Systems U.S.A., Inc.
81 Hartwell Avenue, Suite 300
Lexington, MA 02421 U.S.A.

Global Support

<http://www.fujifilm.com/worldwide/>

Copyright © 1997-2019 FUJIFILM Medical Systems U.S.A., Inc. All rights reserved.

Unauthorized reproduction or distribution of this program, or any portion of it, is strictly prohibited.

Other products and company names mentioned herein may be trademarks of the respective owners.

Certain images/diagrams contained in this manual may not appear clearly on high resolution monitors (1K and 2K)

European Authorized Representative:
FUJIFILM Europe GmbH
Heesenstrasse 31
40549 Düsseldorf, Germany

CE0413

Notified Body:
Intertek Semko AB
Torshamnsgatan 43, Box 1103
SE-164 22 Kista, Sweden



Synapse PACS IHE Integration Statement

Vendor	Product Name	Version	Date
FUJIFILM	Synapse PACS	4.4 and above	09/04/2018
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below.			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Manager / Archive Image Display Performed Procedure Step Manager		
Patient Information Reconciliation	Image Manager / Archive Performed Procedure Step Manager		
Access to Radiology Information	Image Manager / Archive Report Reader		
Simple Image and Numeric Report	Report Reader		
Consistent Time	Time Client		
Teaching File and Clinical Trial Export	Export Selector		
Consistent Presentation of Images	Image Manager / Archive		
Key Image Note	Image Manager / Archive		
Mammography Image	Image Manager / Archive		
Portable Data for Imaging	Portable Media Creator		
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	Basic Patient Privacy Enforcement Basic Patient Privacy Proof	
	Document Source Integrated Document Source/Repository	Basic Patient Privacy Enforcement	
	Document Repository		
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer		
	Imaging Document Source		
	Document Consumer		
	Document Repository		
Audit Trail and Node Authentication	Secure Application		
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer		
Patient Demographics Query	Patient Demographics Consumer		

Basic Patient Privacy Consents	Content Consumer	Basic Patient Privacy Acknowledgement View
Imaging Object Change Management (IOCM)	Change Requester Image Manager / Archive	
Internet address for vendor's IHE information: http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/		
Links to Standards Conformance Statements of the Implementation		
DICOM	http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/	
HL7	http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/	
Links to general information on IHE		
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org

Proprietary Information: For Distribution within FUJIFILM Medical Systems U.S.A., Inc. and its Authorized Affiliates, Distributors, Partners, and Customers Only. Duplication or Distribution of this document without agreement from FUJIFILM Medical Systems U.S.A., Inc. is strictly prohibited.